surefire®

Surefire® Precision Infusion System

>> DYNAMIC. >> FLEXIBLE. >> TRACKABLE.

Instructions for Use	2
Instrucciones de uso.....	4
Instructions d'Utilisation	6
Gebrauchsanleitung.....	8
Istruzioni per l'uso.....	10
Brukerveiledning.....	12
Instruções de utilização	14
Käyttöohje.....	16
Οδηγίες Χρήσης.....	18
Gebbruksaanwijzingen.....	20
Pokyny k použití.....	22
Brugsanvisning.....	24
Instrukcja użycia	26
Инструкции по применению.....	28
Bruksanvisning.....	30
Kullanma Talimatları	32

 **Manufactured by:**
Surefire Medical, Inc.
6272 W. 91st Avenue
Westminster, CO 80031, USA
Customer Service: (888) 321-5212

U.S. Pat. 8,500,775, 8,696,698, 8,696,699 and U.S. and international patents pending.

EMERGO EUROPE
Molenstraat 15
2513 BH The Hague
The Netherlands

DWG-0307-00 R03

Instructions for Use: Surefire® Precision System

Sterile. Sterilized with ethylene oxide gas. Single use only. Do not resterilize.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Patent Pending

Rx
ONLY

CAUTION: CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. FAILURE TO OBSERVE ALL WARNINGS AND PRECAUTIONS MAY RESULT IN COMPLICATIONS.

DESCRIPTION: The Surefire® Precision Infusion System is a 0.021" lumen coaxial microcatheter with the Expandable Tip at the distal end. The Surefire Precision (Figure 1) serves as the conduit for physician-specified agents such as contrast agents, flush solutions, and embolic beads. It is compatible with standard guide wires up to 0.018", and embolic hydrogel particles 500µm or less in size and glass microspheres 110µm or less in size. The Surefire Precision has a PTFE inner liner to provide a lubricious surface for passage of physician-specified agents and other accessory devices. The device is hydrophilically coated. The soft, pliable, funnel-shaped Expandable Tip (Figures 2a and 2b) is sized for use in vessels of 1.5-3.5mm.

FIGURE 1. SCHEMATIC OF THE SUREFIRE PRECISION INFUSION SYSTEM


FIGURE 2A. COLLAPSED TIP



FIGURE 2B. EXPANDED TIP

There are three radiopaque markers located at the distal end of the Surefire Precision to aid in positioning of Expandable Tip. The Expandable Tip can be expanded or collapsed up to 3 times for re-positioning during an interventional procedure by moving the thumb-slide back and forth. When expanded, the Expandable Tip is designed to improve infusion efficiency of compatible embolic agents while maintaining antegrade flow in various size vessels.

The system is provided sterile (EO) and for single patient use.

1. Intended Use: The Surefire® Precision Infusion System is intended for use in angiographic procedures. It delivers radiopaque media and therapeutic agents to selected sites in the peripheral vascular system.

2. Contraindications: The Surefire® Precision Infusion System is not intended for use in the vasculature of the central nervous system (including the neurovasculature) or central circulatory system (including the coronary vasculature).

3. Warnings

- Do not use the product after the "Use By" date specified on the package
- Only physicians possessing sufficient training, skill and experience in the principles, clinical applications, complications, and side effects commonly associated with similar or same techniques should perform vascular interventional procedures.
- Do not use a device where the integrity of the sterile packaging has been compromised.
- Inspect the device prior to use. If the device appears damaged, replace with another device.
- Do not heat or bend or attempt to shape the catheter tip. It may result in abrasion of the hydrophilic coating or damage to the catheter.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
- Do not expand the device close to an implanted vascular device.
- Never collapse, expand or torque the Expandable Tip against resistance as this could cause vessel trauma, device damage or breakage. If resistance is encountered while collapsing/expanding the tip: Flush the infusion lumen and side-arm, and gently attempt to collapse/expand the tip again. If resistance is still encountered, carefully pull the entire infusion system back into the guiding catheter and remove as a unit from the patient.
- Retraction out of the patient with a fully or partially expanded Expandable Tip may result in damage to the Expandable Tip and/or vessel trauma.
- The device is compatible with solutions containing embolic agents, specifically hydrogel beads ≤ 500 µm and glass beads ≤ 110 µm.
- The maximum guide wire diameter for the Surefire Precision is 0.46 mm (0.018").
- Do not use a power injector to infuse agents other than contrast media, as the catheter may become blocked.
- The maximum dynamic injection pressure during power injection should not exceed 8274 KPa / 1200 psi / 82 atm/bar. The static pressure should not exceed 2068 KPa / 300 psi / 20 atm/bar. Static pressure in excess of this maximum may result in device rupture.
- If flow through the catheter becomes restricted, do not attempt to clear the catheter lumen by infusion. Identify and resolve the cause of the blockage or replace the catheter with a new catheter before resuming infusion.

4. Precautions

- Do not store at extreme temperatures and humidity. Avoid direct sunlight.
- This device is intended for one-time use only. Do not re-sterilize and/or reuse.
- Handle the device carefully to avoid contact with sharp instruments, plastic fittings, or abrasive surfaces that may damage the materials.
- Maintain a continuous heparinized saline flush in order to achieve optimal device performance and prevent or reduce the risk of thrombus formation on the catheter and expandable tip. The use of systemic heparinization to reduce the risk of thromboembolic complications should be considered for some patients, such as those with known hypercoagulable conditions.
- The surface of the distal portion of the microcatheter and the Expandable Tip must be completely wet with heparinized saline to remain lubricious. Keep these areas wet prior to and during use. Minimize handling of the expandable tip to avoid damaging it.
- Do not expose the delivery system to organic solvents (e.g. alcohol) as structural integrity and / or function of the device may be impaired.

• There are three radiopaque markers located at the distal end of the Surefire Precision to aid in positioning the Expandable Tip.

- Do not attempt to advance or torque the Surefire Precision while the Expandable Tip is expanded as this may result in damage to the Expandable Tip and/or vessel trauma.

5. Recommended Accessories

- Appropriate guiding catheter (with inner diameter of 0.046" minimum) and compatible sheath introducer to accommodate the Surefire Precision.
- 0.014" (0.36 mm), 0.016" (0.41 mm), 0.018" (0.46mm) straight tip or curved tip, guide wire
- Heparinized saline or equivalent flushing solution
- 1cc, 3cc, and 10cc or 20cc luer lock syringes

6. Adverse Events: Complications specific to Surefire Precision procedures may occur at any time during or after the procedure and may include, but are not limited to, the following:

- Access site complications
- Distal embolization of system components
- Vessel dissection, perforation, rupture, and hemorrhage
- Vascular Thrombosis
- Embolism
- Ischemia
- Infection
- Allergic reaction
- Vasospasm

7. Preparations for Use

- Remove the Surefire Precision sterile pouch from the device box and inspect for damage. If there is any damage to the pouch or the device, discard and replace with a new sterile device.
- Place the appropriate sheath introducer and guiding catheter using standard percutaneous technique. A Tuohy-Borst adapter should be attached to the guiding catheter to allow for continuous flush of the guiding catheter with heparinized saline.

Surefire Precision Preparation:

- Carefully open the sterile pouch. Employ aseptic technique during removal of the device from the packaging and during use.
- Gently remove the card with the coiled dispenser tube containing the Surefire Precision.
- Fill a 1cc, 3cc and a 10cc or 20cc syringe (as indicated on the packaging card) with heparinized saline. Prepare device as follows:
 - 1cc: Flush the side arm. Confirm that liquid exits at the distal section. Refill and keep the luer lock syringe attached to the sideport.
 - 3cc: Flush the infusion port.
 - 10cc or 20cc (as indicated on the packaging card): Flush the distal end of the coiled dispenser tube in order to hydrate the infusion system.
- Move the orange lever to the unlock position. Advance the thumb-slide to fully collapse the expandable tip (see Figures 3a and 3b).

FIGURE 3A. LEVER IN UNLOCK POSITION




FIGURE 3B. THUMB-SLIDE

- Move the orange lever to the lock position (see Figure 4)

FIGURE 4. LEVER IN LOCK POSITION

- Carefully remove the catheter from the dispenser tube, taking care not to damage the catheter upon removal.
- Remove the orange tip protector.
- Inspect the device and expandable tip thoroughly to ensure that it is not kinked or otherwise damaged. If there is any damage, replace with a new sterile device.
- Prepare an appropriately sized guide wire according to the manufacturer's instructions for use. Insert the guide wire into the Surefire Precision.

8. Directions for Use

Positioning

- Insert the Surefire Precision with the guide wire into the guiding catheter. Tighten the Tuohy-Borst adapter around the microcatheter to prevent backflow while still allowing for movement of the microcatheter through the Tuohy-Borst, avoid over-tightening.
 - Note:** Handle the distal end of the Surefire Precision with care. Distortion or flattening of the device may prevent the Surefire Precision from advancing into the lumen of the guiding catheter.
- With Tip in collapsed position, track the Surefire Precision to the target location.
 - Note:** There are three radiopaque markers on the distal end of the Surefire Precision. See Figure 2.
 - Caution:** Do not advance or torque the system against resistance. If the cause of resistance cannot be determined, carefully withdraw catheter.
 - Note:** It is important to adequately flush the infusion lumen throughout the procedure to prevent the backflow of blood into the device lumen.

Tip Expansion and Diagnostic / Embolic Infusion

- Move the orange lever to the unlock position.
- Flush the sideport before and after each opening and closing of the tip.
- Use fluoroscopic guidance when opening the tip, move the thumb-slide slowly to expand the tip. The tip is fully deployed when all three marker bands are equally spaced.
 - Note:** It may not be necessary to pull the thumb-slide to its proximal limit in order to fully deploy the expandable tip. Fluoroscopic imaging of the markers should be used to determine tip position.
- Move the orange lever on the handle to the lock position.
- Remove the guide wire. Inject with contrast to confirm infusion location and apposition of tip to vessel wall under fluoroscopy. If repositioning is needed, flush the sideport and collapse the tip prior to repositioning.
 - Note:** Always inject contrast through the infusion port to confirm location prior to infusion. Confirm that contrast can be injected freely without excess pressure prior to infusion of embolic agents.
- Once position in the desired location is confirmed, close the Tuohy-Borst on the guiding catheter to prevent movement of the Surefire Precision. Do not overtighten.
 - Caution:** Do not attempt to advance or torque the Surefire Precision while the Expandable Tip is expanded as this may result in damage to the Expandable Tip and/or vessel trauma.
- Confirm antegrade flow. If antegrade flow is less than desired, flush sidearm and infusion hub of device, reposition as needed and reassess antegrade flow.
- Adequately flush the device with heparinized saline prior to initial introduction of therapeutic agents.
- Infuse diagnostic, embolic or therapeutic agents according to the manufacturer's instructions for Use.

Device Retraction

- Flush the Tuohy-Borst on the guiding catheter.
- Flush the Surefire sideport.
- Open the Tuohy-Borst on the guiding catheter.
- Move the orange lever on the Surefire handle to the unlock position.
- Under fluoroscopic guidance, advance the inner microcatheter forward to collapse the tip (by advancing the thumb-slide) while simultaneously retracting the system as a whole. Lock the handle, flush the sideport, remove.
 - Warning:** Retraction of a fully deployed Expandable Tip may result in damage to the Expandable Tip and/or vessel trauma.
 - Warning:** Never retract the Expandable Tip against resistance as this could cause vessel trauma and/or device damage or breakage. If resistance is encountered when collapsing the tip, stop collapsing the tip and carefully pull the entire infusion system back into the guiding catheter and remove as a unit from the patient.
- If the Precision will be used again within the same procedure, flush the infusion port and sidearm, and keep the device soaked in heparinized solution.
- After use, dispose in accordance with hospital, administrative and/or local government policy.
 - Caution:** This device is intended for single patient use only. Do not re-sterilize and/or reuse.

Instructions for Using a Power Injector

A power injector can be used to infuse contrast media through the device. The flow rate depends upon factors such as the viscosity of the contrast media, which varies with the type and temperature of the media, the model and setting of the power injector, and how the injector is connected to the device. The observed flow rate values indicated below are for reference only.

Catheter ID	Usable Length	Infusion Medium	Viscosity (cP)	Set Condition		Actual Flow Rate (mL/sec)			Dead Space Volume (mL)
				Flow Rate (mL/sec)	Volume (mL)	4137 kPa (600 psi)	6205 kPa (900 psi)	8274 kPa (1200 psi)	
0.021"	120 cm	0.9% Saline	1.0	3.0	10.0	2.0	3.0	3.0	0.37
		100% Contrast (iodine 370 mg/mL)	10.5	3.0	10.0	0.9	1.4	1.9	
	150 cm	0.9% Saline	1.0	3.0	10.0	1.7	2.8	3.0	0.44
		100% Contrast (iodine 370 mg/mL)	10.5	3.0	10.0	0.6	1.0	1.4	

Injector Used: MEDRAD Mark V Provis Contrast Media Temperature: 37 °C

Label Symbols

Heparinized saline flush volume indicated.
 Thumbslide unlocked on left.
 Thumbslide locked on right.
 Do not use if package is damaged.
 Does not contain natural latex rubber.
 Recommended Vessel Size
 Inner Diameter
 Maximum Guidewire
 Minimum Guide Catheter Inner Diameter
 Maximum Dynamic Injection Pressure

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY
 THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE SUREFIRE MEDICAL, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL SUREFIRE MEDICAL, INC. BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND SUREFIRE MEDICAL, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN DESCRIPTIONS OR SPECIFICATIONS IN SUREFIRE, INC PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES. SUREFIRE MEDICAL, INC. WILL NOT BE RESPONSIBLE FOR ANY DIRECT OR CONSEQUENTIAL DAMAGES RESULTING FROM REUSE OF THE PRODUCT.

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Specification and Inspection Note Table		
Specification	Inspection Method	Insp. Req'd.
1. Packaging		
A. Components will be packaged in a way not to compromise the integrity of the components.	Visual inspection (unaided eye)	Yes
2. Material Specification(s)		
A. Trim Size: 18 x 9 inches; folds/stitches to 9 x 9 inches B. No. of Pages: Pages numbered 2 through 33, plus self-cover C. Paper: 40# White Opaque Offset or 50# White Opaque Digital D. Printing Color: black, no bleeds E. Finishing: Fold, stitch and trim to 9 x 9 inches. Carton packed.	Visual inspection (unaided eye) (AQL=4.0)	Yes
3. Identification and Traceability		
A. Packing list indicating (at a minimum) the PO #, Surefire Medical's P/N and revision, description and quantity in each package.	Attached to or inside shipping container	Yes
B. Bulk and Assembly Packaging will be labeled with (at a minimum) Surefire Medical's P/N and revision, description, and quantity in each package.	Visual inspection (unaided eye)	Yes
4. Physical / Functional Specifications		
A. Free from obvious discoloration, contamination, and other damage.	Visual inspection (unaided eye) (AQL=4.0)	Yes
B. Inspect one random piece from vendor lot to ensure that it matches the component drawing.	Visual inspection (unaided eye)	Yes
5. Dimensional Inspection		
N/A	N/A	N/A
6. Supplier Processing Specifications		
N/A	N/A	N/A