# surefire®

## Surefire® Flexion Guiding Catheter

Instructions for Use .....	2
Instrucciones de uso.....	4
Instructions d'utilisation.....	6
Gebrauchsanleitung.....	8
Istruzioni per l'uso.....	10
Instruksene for bruk.....	12
Instruções de utilização.....	14
Käyttöohjeet.....	16
Οδηγίες Χρήσης.....	18
Instructies voor Gebruik.....	20
Pokyny k použití.....	22
Brugsanvisning.....	24
Instrukcja użycia.....	26
Инструкции по применению.....	28
Bruksanvisning.....	30
Kullanma Talima.....	32

surefire®  
MEDICAL INC.

Surefire® Medical, Inc.  
6272 W. 91st Ave.  
Westminster, CO 80031 U.S.A.  
Customer Service: (888) 321-5212

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

# Surefire® Flexion Guiding Catheter Instructions for Use

ENGLISH

**Sterile. Sterilized with ethylene oxide gas. Single use only. Do not resterilize.**

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Patent Pending**

**CAUTION: CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. FAILURE TO OBSERVE ALL WARNINGS AND PRECAUTIONS MAY RESULT IN COMPLICATIONS.**

**DESCRIPTION:** The Surefire® Flexion Guiding Catheter provides a pathway to introduce and facilitate the advancement of devices into the peripheral vascular system. The Surefire® Flexion Guiding Catheter is a single-lumen, braided catheter with a soft distal tip and a proximal Luer-Lock hub. The Surefire® Flexion Guiding Catheter is available in a variety of pre-shaped tip designs to accommodate access and positioning in a range of peripheral vascular anatomies.

**1. Intended Use:** The Surefire® Flexion Guiding Catheter is intended to provide a pathway through which therapeutic devices are introduced. The Surefire® Flexion Guiding Catheter is intended to be used in the peripheral vascular system.

**2. Contraindications:** The Surefire® Flexion Guiding Catheter is not intended for use in the vasculature of the central nervous system (including the neurovasculature) and central circulatory system (including the coronary vasculature).

### 3. Warnings

- Do not use the product after the "Use By" date specified on the package.
- Only physicians possessing sufficient training, skill and experience in the principles, clinical applications, complications, and side effects commonly associated with similar or same techniques should perform vascular interventional procedures.
- Do not use a device where the integrity of the sterile packaging has been compromised.
- Inspect the device prior to use. If the device appears damaged, replace with another device.
- Do not heat or bend the catheter tip. It may result in damage to the catheter.
- Do not position the device close to an implanted vascular device.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter if resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Never advance, retract, or torque the device in the vessel against resistance as this could cause vessel trauma or device damage.
- Always use a guide wire to remove the catheter from the vasculature. Failure to do so may result in damage to the vessel, puncture site, product, or all three.
- The maximum guide wire diameter for use with the catheter is 0.038" (0.97 mm).
- During power injection, do not exceed the maximum dynamic injection pressure of 600 psi (4137 kPa). The static pressure should not exceed 300 psi (2068 kPa). Static pressure in excess of this maximum may result in device rupture.
- Flow rates have been established using OXILAN® 300 (9.4 cPs @ 20°C) for a 2 second injection. Since a number of variables (such as media viscosity and temperature) can produce different flow rates, this information should be used for reference only.

### Generated Pressures at Given Flow Rates

ml/sec	4	6	8	10	15
psi	111	190	278	369	599
kg/cm2	8	13	20	26	42

### 4. Precautions

- Do not store at extreme temperatures and humidity. Avoid direct sunlight.
- This device is intended for one-time use only. Do not re-sterilize and/or reuse.
- Handle the device carefully to avoid contact with sharp instruments, plastic fittings, or abrasive surfaces that may damage the materials.
- Do not torque the catheter if it is kinked. Doing so may cause damage to the shaft that could result in catheter separation. If the catheter shaft becomes kinked, withdraw the guiding catheter, guide wire, and catheter sheath introducer.
- Consider the use of systemic heparinization to prevent or reduce the possibility of thrombus formation on the surface of the catheter.
- Do not expose the catheter to organic solvents as structural integrity and/or function of the device may be impaired.
- The distal tip of the Catheter is comprised of a radiopaque polymer that should be visible for accurate placement. The proximal shaft segment is also radiopaque.

### 5. Recommended Accessories

- 5F x 10 cm sheath introducer with side flushing port
- 0.038" (0.97 mm) x 150 cm, guide wire
- Heparinized saline or equivalent flushing solution
- Tuohy-Borst

**6. Adverse Events:** Complications specific to the Surefire® Flexion Guiding Catheter procedures may occur at any time during or after the procedure and may include, but are not limited to, the following:

- Access site complications
- Vessel dissection, perforation, rupture, and hemorrhage
- Vascular Thrombosis
- Embolism
- Ischemia
- Infection
- Allergic reaction
- Vasospasm

## 7. Preparations for Use

### Access Site and Treatment Site Selection:

- Access the target vessel using standard percutaneous technique.
- Insert the recommended sheath introducer using standard technique.
- Secure the sheath introducer in place.

### Catheter Preparation:

- Remove the catheter sterile pouch from the device box and inspect for damage. If there is any damage to the pouch, discard and replace with new sterile device.
- Employ aseptic techniques during removal from the packaging and during use.
- Carefully open the catheter sterile pouch and gently remove the catheter from the dispenser card. Inspect thoroughly, ensuring that it is not kinked or otherwise damaged. If there is any damage, discard and replace with a new device.
- Flush the catheter with heparinized saline, and then attach a Tuohy-Borst to the hub of the catheter.

## 8. Surefire® Flexion Guiding Catheter Directions for Use

### Catheter Introduction:

- Place a catheter sheath introducer using standard percutaneous technique.
- Insert the Surefire® Flexion Guiding Catheter into the vasculature through the sheath introducer and over an appropriately sized guide wire.
- Place the Catheter over the guide wire and advance the catheter in the normal manner.
- Manipulate the Catheter carefully into the desired artery using fluoroscopy.
- Remove the guide wire from the Catheter prior to the introduction of other devices or contrast agents.
- Flush the Catheter with heparinized saline periodically, and prior to the introduction of other devices.

**Warning:** Do not advance or torque the system against resistance. If the cause of resistance cannot be determined, withdraw catheter and guide wire together as a unit to prevent potential vessel and product damage.


**Note:** The Catheter shaft and distal tip contain radiopaque polymers that should be visible for accurate placement.

### Catheter Retraction:

- Under fluoroscopic guidance, retract the Catheter.  
**Warning:** Never advance, retract, or torque the Catheter against resistance as this could cause device damage.
- After use, dispose in accordance with hospital, administrative and/or local government policy.

**Warning:** This device is intended for one-time use only. Do not re-sterilize and/or reuse.

## Label Symbols

## DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE SUREFIRE® MEDICAL, INC. PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL SUREFIRE® MEDICAL, INC. BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND SUREFIRE® MEDICAL, INC. TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

DESCRIPTIONS OR SPECIFICATIONS IN SUREFIRE®, INC PRINTED MATTER, INCLUDING THIS PUBLICATION, ARE MEANT SOLELY TO GENERALLY DESCRIBE THE PRODUCT AT THE TIME OF MANUFACTURE AND DO NOT CONSTITUTE ANY EXPRESS WARRANTIES.

SUREFIRE® MEDICAL, INC. WILL NOT BE RESPONSIBLE FOR ANY DIRECT OR CONSEQUENTIAL DAMAGES RESULTING FROM REUSE OF THE PRODUCT.

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6272 W. 91st Avenue • Westminster, CO 80031 • U.S.A.